

**Product Name :**  
Newborn Intubation Baby Model

**Product Code :**  
NLAB-MEDICALMODELAB11024



**Description :**

Newborn Intubation Baby Model

**Technical Specification :**

The size and configuration of this model is almost as same as the normal newborn baby.

The head, upper limbs, lower limbs and all joints are flexible.

It can demonstrate the characteristic of the newborn baby and

all kinds of operations in clinic education: includes different

parturition method, newborn baby bath etc.

Made of soft PVC plastic.

Washable.

Length: 52cm, head diameter: 34 cm.

This model is used to demonstrate the newborn baby.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003