

Product Name :
Obd 2 Code Reader - OBD II

Product Code :
TK111



Description :

Obd 2 Code Reader - OBD II

Technical Specification :

This compact hand-held OBD2 tool offers light do-it-yourselfers and consumers a quick & easy way to find out why the "CHECK ENGINE" light is illuminated on 1996 and newer cars, light trucks, mini-vans and SUVs - foreign and domestic.

This little tool packs a big punch with its patented single-screen display and system status LEDs for a quick check.

Large easy-to-read LCD screen.

Retrieves both Generic (P0) and Manufacturer specific (P1) diagnostic codes.

Can LINK and retrieve codes in less than 10 seconds.

On a single digital screen, the tester provides information of:

Monitor Status / I/M Readiness to determine if vehicle is ready for an emissions test.

MIL (Malfunction Indicator Lamp) status -- (1st code in sequence shows which code set the "Check Engine" light to go on)

Current Diagnostic Trouble Codes (DTCs)

Alerts of any Pending Codes

Auto-Monitoring function periodically retrieves the latest code and I/M readiness status from a vehicle's computer. This gives the user the ability to verify if the repair was completed properly and see if all the Monitors completed their diagnostic testing.

ERASE button clears stored computer codes, Monitor status and turns off "Check Engine" warning lamp.

Includes a 6-foot removable and serviceable Data Link Connector (DLC) cable.

Includes a detailed instruction manual in English, French and Spanish with in-depth information on OBD II

systems and codes.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003