

Product Name : Ohm Meter Measure Resistance	Product Code : NLAB-TESTINGLABM37003
Description : Ohm Meter Measure Resistance	
Technical Specification : Educational Desk Stand Ohm's meter Ohm's meter is designed for measuring the Resistance in electric circuits It is housed in unbreakable ABS case with clear transparent front cover and operated on 1.5V AA cell. Size : 100 x 100 x 140mm	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003