

**Product Name :**  
Oil Insulation Tester Machines

**Product Code :**  
NLAB-TESTINGLABM23001



**Description :**

Oil Insulation Tester Machines

**Technical Specification :**

Fully automatic measurement of the dielectric strength of liquid insulations up to 100 kV

Excellent design for high efficiency and usability

Extraordinarily powerful for continuous use in the laboratory

Clear and reliable results through the latest measurement technology

Repeatable breakdown measurement in mineral oil, silicone and ester fluids

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003