

Product Name :
Ointment Cream Filling Machine

Product Code :
NLAB-SCIENTIFICEQ28014



Description :

Ointment Cream Filling Machine

Technical Specification :

Durable

Weather proof

Temperature resistance

Technical Specifications:

Power Supply: 220/110 V 50/60 Hz

Filling Speed: 5-20 Bottles/min/head

Filling Accuracy: $\leq \pm 1\%$

Air Pressure: 0.4-0.6 MPa

Filling Range: 25-60 ml 50-125 ml 100-250 ml 200-500 ml

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003