

**Product Name :**  
Open Bedside Locker

**Product Code :**  
NLABMEDICAL-M2500003



**Description :**

Open Bedside Locker

**Technical Specification :**

Specification :-

Frame Work of Open Bedside Locker is made of Square tube.

It has Stainless Steel tray top.

Open Bedside Locker is Pre Treated Epoxy Powder Coated.

Product Size : 36 L x 36 W x 81 H cms.

We are Manufacturer, Exporters and Suppliers of Open Bedside Locker at tender and export competitive prices from India.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

