

Product Name : Operation Theatre Light	Product Code : NLABMEDICAL-M2750002
Description :	
Operation Theatre Light	
Technical Specification :	
Dome Head: 521mm	
Halogen Bulbs: 24V x 150 W x 2 nos	
Light Intensity at 1 mt: 50,000 Lux	
Dichoric Glass Reflector: 2000 mm	
Color Temperature: 4200 K \pm 300	
Temp rise in field: 2°-3°C	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003