

Product Name :
Opposed-piston Engine Cutaway

Product Code :
NLAB-AUTOMOTIVETVET19003



Description :

Opposed-piston Engine Cutaway

Technical Specification :

This cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts, cross-sections, lubricating circuits, fuel system, cooling system etc. Many parts have been chromium-plated and galvanized for a longer life.

Main Technical Specifications:

- Lycoming/Piaggio/Continental 6-opposed cylinders engine
- Air cooling
- Gear distribution with camshaft in the crankcase
- Ignition with magneto
- Single-body carburettor

Automotive Technologies, Cutaways Models, Aeronautics Engines Cutaways Equipments for Automotive Training Lab

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003