

**Product Name :**  
Optical Bench Candle Holder

**Product Code :**  
NLAB-EDUSCHOOLEQ174097



**Description :**

Optical Bench Candle Holder

**Technical Specification :**

Cylindrical 20 mm internal diameter x 22 mm deep painted matt black is mounted coaxially upon a polished brass rod 35 x 6 mm diameter.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003