


Product Name :
Optical Bench Set

Product Code :
NLAB-SCHOOLABS19003



Description :

Optical Bench Set

Technical Specification :

Study the image formation and the focal length of the resulting image.

The equipment within the set includes Scale rail kit The length is about 1 meter and can be moved easily. Per use, convex-flat lens($\phi=50\text{mm}$, $F=300\pm 2\text{mm}$), convex lens($\phi=40\text{mm}$, $F=100\pm 2\text{mm}$),

convex lens($\phi=30\text{ mm}$, $F=50\pm 2\text{mm}$), concave lens ($\phi=30\text{mm}$, $F=-75\pm 4\text{mm}$), glass insert – opaque

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003