


**Product Name :**  
Optical Bench Type 2

**Product Code :**  
NLAB-SCHOOLABS19004



**Description :**

Optical Bench Type 2

**Technical Specification :**

Study the image and the focal length. It can be used to calculate image size and direction of image formation. The kit includes: light source 220 V 40 W, various sizes of lenses, object plate, bracket, metal rail with base and instruction manual.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003