


Product Name :
Optical Bench Type 2

Product Code :
NLAB-SCHOOLABS19004



Description :

Optical Bench Type 2

Technical Specification :

Study the image and the focal length. It can be used to calculate image size and direction of image formation. The kit includes: light source 220 V 40 W, various sizes of lenses, object plate, bracket, metal rail with base and instruction manual.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003