

Product Name :
Orested Apparatus

Product Code :
NLAB-EDUSCHOOLEQ174098



Description :

Orested Apparatus

Technical Specification :

Rested's Apparatus for Physics Lab To show the action of electric current on a magnet. Magnetic needle 75 mm long with brass bearing supported by pivot mounted on a plastic base & surrounded by a rectangular copper rod fitted with terminals. Dimensions 80 x 40 mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003