

Product Name :
Nauga Lab Organ Bath

Product Code :
NLAB-LABORATORYEM21014



Description :

Nauga Lab Organ Bath

Technical Specification :

Made of transparent acrylic sheet, chemically welded to make it perfectly leak-proof.

Capacity of bath 750ml.

Complete with uprights, organ & oxygen tubes, capillary tubing, frontal lever.

Temperature by immersion heater controlled by a Thermostat.

Chamber Size: 8" x 5" x 5"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003