

Product Name :
Organ Jar Bath Digital

Product Code :
NLAB-EDUSCHOOLEQ46004



Description :

Organ Jar Bath Digital

Technical Specification :

Made of transparent acrylic sheet, chemically welded to make it perfectly leak-proof.

Capacity of bath 750 ml.

Complete with uprights, organ & oxygen tubes, capillary tubing, frontal lever.

Temperature by immersion heater controlled by a microprocessor based digital temperature controller

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003