

<b>Product Name :</b> Naugra Lab Organ Jar Baths Digital	<b>Product Code :</b> NLAB-LABORATORYEM21015
<b>Description :</b>	
Naugra Lab Organ Jar Baths Digital	
<b>Technical Specification :</b>	
Made of transparent acrylic sheet, chemically welded to make it perfectly leak-proof.	
Capacity of bath 750 ml.	
Complete with uprights, organ & oxygen tubes, capillary tubing, frontal lever.	
Temperature by immersion heater controlled by a microprocessor based digital temperature controller	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003