

**Product Name :**  
Osborne Reynolds Flow Regimes Demonstration Module

**Product Code :**  
NLAB-EDUCATIONTVET14006



**Description :**

Osborne Reynolds Flow Regimes Demonstration Module

**Technical Specification :**

The set makes it possible to carry out the classic Osborne Reynolds experiment and to observe laminar, transient, turbulent flows and the velocity profile.

The test tube has a diameter of 10 mm and a length of 700 mm. The flow through this tube is regulated by a valve. The injection system

of dye is similar to that used by Reynolds (0.45 l capacity tank)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003