

Product Name :
Oscilloscope

Product Code :
ES848



Description :

Oscilloscope

Technical Specification :

Digital 25MHz.
Sampling frequency, 100MS/s f=50Hz.
Voltage =100-240V
Reading-out with the cursor.
Five automatic measurement functions.
Color liquid crystal display of high resolution.
Storage and call-out of waveforms.
Automatic setting function provided capable of fast setting.
Multiple-waveform calculation function.
Implementation of detecting the average and peak values of the waveform.
Digital real-time oscilloscope.
Edge and video triggering function.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

