

**Product Name :**  
Overbed Bedside Table Adjustable 2 Section

**Product Code :**  
NLABMEDICAL-M2810021



**Description :**

Overbed Bedside Table Adjustable 2 Section

**Technical Specification :**

Frame is made of Rectangle & Square Tube.

The body section is mounted on 5cms castors.

The size of the Laminated top is 80L x 40W cms.

The Laminated top comes in two sections.

The height of the Overbed Bedside table is adjustable manually by Gear Handle.

Knock Down Construction to save freight.

We are Manufacturer, Exporters, and Suppliers of Overbed Bedside Table Adjustable 2 Section at tender and export competitive prices from India.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments