

Product Name :
Overbed Bedside Table Adjustable 2 Section

Product Code :
NLABMEDICAL-M2690004



Description :

Overbed Bedside Table Adjustable 2 Section

Technical Specification :

Specification :-

Frame is made of Rectangle & Square Tube.

The body section is mounted on 5cms castors.

The size of Laminated top is 80L x 40W cms.

The Laminated top comes in two sections.

The height of Overbed Bedside table is adjustable manually by Gear Handle.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Overbed Bedside Table Adjustable 2 Section at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003