

Product Name : Oxygen Concentrator	Product Code : NLABMEDICAL-M3310003
Description :	
Oxygen Concentrator	
Technical Specification :	
Offers a broad therapy spectrum with 0.5 to 5 liters O2/min.	
Mobile dose monitor with its 20 meters long hose permits greater freedom of movement.	
Very quite -only 40 dB.	
Oxygen output more than 90%, outlet pressure 0.6 bar.	
Alarms for mains failure/device failure and low O2 concentration.	
Small, light weight (only 20 kgs.), modular construction, snap closures, simplest filter change & low maintenance.	
High performance and output even at 2000 meters above.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003