

Product Name :
Pal and Gal Programmer Trainer

Product Code :
NLAB-ENGINEERINGLB54042



Description :

Pal and Gal Programmer Trainer

Technical Specification :

It features a universal intelligent programmer, PC computer, interface, and software with numerous drivers for many brands of PLDS and also a text book with exercises.

This trainer comprises a complete system for the training of Programmable Logic Devices (PLD).

The Programmer is a universal device programmer, which works through the PC's parallel port.

Assortments of PLDS are also provided for experimentation.

We are well-known manufacturers, OEM suppliers of Pal and Gal Programmer Trainer for Applied Machines. Contact us for high quality Pal and Gal Programmer Trainer for Applied Machines for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003