

Product Name :
Paraffin Dispenser

Product Code :
NLAB-EDUSCHOOLEQ29002



Description :

Paraffin Dispenser

Technical Specification :

Paraffin Dispenser

- Embedding bath-single type dispenser which is double walled, both made of stainless steel sheet, gap filled with glass wool for temp. Insulation.
- Paraffin Dispensers have high heating performance and provide precise Temperature range ambient to **70° C** controlled by hydraulic type thermostat.
- Completely supplied with Pilot Lamp, Temp. Control Knob, Cord-Plug etc.
- To work on **220V AC 50Hz**, single phase.

Optional:

- Stainless steel basket.

Available in capacity:

1 litre

2 litre

5 litre

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003