

**Product Name :**  
Paraffin Dispenser

**Product Code :**  
NLAB-EDUSCHOOLEQ29002



**Description :**

Paraffin Dispenser

**Technical Specification :**

**Paraffin Dispenser**

- Embedding bath-single type dispenser which is double walled, both made of stainless steel sheet, gap filled with glass wool for temp. Insulation.
- Paraffin Dispensers have high heating performance and provide precise Temperature range ambient to **70° C** controlled by hydraulic type thermostat.
- Completely supplied with Pilot Lamp, Temp. Control Knob, Cord-Plug etc.
- To work on **220V AC 50Hz**, single phase.

**Optional:**

- Stainless steel basket.

**Available in capacity:**

1 litre

2 litre

5 litre

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003