

Product Name :
Paraffin Embedding Bath

Product Code :
NLAB-EDUSCHOOLEQ29012



Description :

Paraffin Embedding Bath

Technical Specification :

Paraffin Embedding Bath

The body is made of mild steel sheet which is duly painted.

The top is made of M.S. / S.S. sheet.

Completed with wax pots with covers.

Temperature ambient is upto 70° C controlled by a thermostat to work on 220V AC 50Hz single phase.

The following compositions are available:-

6 GURSS TUBES

8 GURSS TUBES

8 GURSS TUBES

12 GURSS TUBES

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003