

Product Name :
Particle Size Analysis Set Hydrometer Method

Product Code :
NLAB-SOILELAB140028



Description :

Particle Size Analysis Set Hydrometer Method

Technical Specification :

This equipment is used to determine the quantitative size distribution of very fine particle in soils such as clay and silt.

The complete set comprises:

Sedimentation Cylinder, 1000 ml capacity (6pieces) with rubber bungs. Soil Hydrometer, 151H, range 0.995 to 1.038g/ml with div. 0.001. Glass Tank, dimensions: 600x300x380 mm

Heater complete with thermostat, circulation unit. 230V Single Phase 50Hz 1000W

Glass Thermometer, range -10 to 50° C., Sub-div. 0.5° C.

Beaker, Borosilicate Glass, 250 ml capacity

Sodium Hexametaphosphate, 1000 g

High Speed Stirrer, 4000 R.P.M, complete with baffle cup, for dispersing soilparticles in water

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003