

Product Name :
Particle Size Analysis Test Set

Product Code :
NLAB-VOCATIONALTVET18009



Description :

Particle Size Analysis Test Set

Technical Specification :

Capacity: 8 hydrometer jars.

Construction: Stainless steel tank; enamel finish exterior.

Jar Support: Stainless steel grid.

Temperature Range: Ambient to 60°C (140°F).

Sensitivity: $\pm 0.2^{\circ}\text{C}$ at 20°C (68°F).

Circulating System: Drive pump with motor.

Heating System: Immersion type; 700-watt heater in stainless steel mixing chamber.

Plumbing Ports: Overflow; Drain; Inflow and outflow for the cooling system.

Inside Dimensions: 959 x 197 x 419 mm (37-3/4" w. x 7-3/4" d. x 16-1/2" h.).

Outside Dimensions: 1,016 x 254 x 584 mm. (40" w. x 10" d. x 23" h.).

Weight: Net 36 kg 80 lbs.

110vAC to 230vAC, 50/60 Hz.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003