

Product Name : Pascal's Ball Metallic	Product Code : NLAB-EDUSCHOOLEQ1730168
Description :	
Pascal's Ball Metallic	
Technical Specification :	
To demonstrate the equal transmission by liquids of pressure in all direction to prove pascal's law. The ball is joined to piston and cylinder. Comprising a hollow metal with orifices. Size: 33 x 11 x 11	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003