

<b>Product Name :</b> Pascal's Ball Metallic	<b>Product Code :</b> NLAB-EDUSCHOOLEQ1730168
<b>Description :</b> Pascal's Ball Metallic	
<b>Technical Specification :</b> To demonstrate the equal transmission by liquids of pressure in all direction to prove pascal's law.  The ball is joined to piston and cylinder.  Comprising a hollow metal with orifices.  Size: 33 x 11 x 11	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003