

**Product Name :**  
Pascals Vases Apparatus

**Product Code :**  
NLAB-EDUSCHOOLEQ174099



**Description :**

Pascals Vases Apparatus

**Technical Specification :**

For use as a light operated timing gate. It consists of a photo detector and light source, housed in two plastic cases. 4 mm sockets are provided for connecting it to electronic timer gate.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003