

Product Name :
Pascals Vases Apparatus

Product Code :
NLAB-EDUSCHOOLEQ174099



Description :

Pascals Vases Apparatus

Technical Specification :

For use as a light operated timing gate. It consists of a photo detector and light source, housed in two plastic cases. 4 mm sockets are provided for connecting it to electronic timer gate.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003