


Product Name :
Pat Test Brass Ring Mould

Product Code :
NLAB-CEMETLAB130011



Description :

Pat Test Brass Ring Mould

Technical Specification :

The mould is 100 mm dia. x 5 mm deep with an internal taper of 5°. It is supplied complete with glass plate. Three are required for each test.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003