

**Product Name :**  
PCB Curing Oven Machine

**Product Code :**  
NLAB-ELECTRONICSAB180002



**Description :**

PCB Curing Oven Machine

**Technical Specification :**

PCB Curing Oven Machine For Electrical Lab Training.

PCB Curing Oven Machine  
Tabletop unit for curing of liquid protoresist or drying of PCB's.  
Timer controlled system allows fast and efficient PCB curing.  
PCB size 10" x 12"  
Timers : 0-30 minutes.  
Temperature : 0-90 degree C

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003