

Product Name :
PCB Shearing Machine

Product Code :
NLAB-ELECTRONICSAB180004



Description :

PCB Shearing Machine

Technical Specification :

PCB Shearing Machine For Electrical Lab Training.

PCB Shearing Machine

A Bench top shearing machine to PCB's and Laminators.

Structure with central handle allows ease of cutting.

Width : 300mm.

Table size : 400mm x 350mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003