

**Product Name :**  
PCB Shearing Machine

**Product Code :**  
NLAB-ELECTRONICSAB180004



**Description :**

PCB Shearing Machine

**Technical Specification :**

PCB Shearing Machine For Electrical Lab Training.

PCB Shearing Machine

A Bench top shearing machine to PCB's and Laminators.

Structure with central handle allows ease of cutting.

Width : 300mm.

Table size : 400mm x 350mm

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003