

Product Name :
Pediatric Bed

Product Code :
NLABMEDICAL-M2800003

Description :

Pediatric Bed India, Pediatric Bed manufacturer India, suppliers, Medical Pediatric Bed India, Hospital Pediatric Bed India

Technical Specification :

Pediatric Bed is very useful for patients who met with a severe or major accident

This is specially designed for the injured patients who cannot lie down straight on the bed

This bed is furnished with hard and precise steel rod to give an additional support to the patient

These beds are being used in almost in every hospital

The height of the upper round rod is connected equally to give support to its lower stumps

Frame work made of rectangular M.S. tubes

Perforated CRC Sheet duly machine press bent

M.S. Powder Coated Tubular Head & Foot Bows of equal height

Drop down Safety side railings on both sides

Bed Mounted on PVC Stumps

Pre Treated and Epoxy powder Coated

Frame work of precise steel tubes

Top of perforated CRCA sheet main frame of rectangular pipe with full length drop side railing Bows of equal height mounted on rubber stump Finish in epoxy powder coated

approximate Size{LxWxH}

{137x 76 x60} mm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003