

| | |
|---|---|
| Product Name : Pediatric Bed | Product Code : NLABMEDICAL-M2800003 |
| Description : Pediatric Bed India, Pediatric Bed manufacturer India, suppliers, Medical Pediatric Bed India, Hospital Pediatric Bed India | |
| Technical Specification : Pediatric Bed is very useful for patients who met with a severe or major accident This is specially designed for the injured patients who cannot lie down straight on the bed This bed is furnished with hard and precise steel rod to give an additional support to the patient These beds are being used in almost in every hospital The height of the upper round rod is connected equally to give support to its lower stumps Frame work made of rectangular M.S. tubes Perforated CRC Sheet duly machine press bent M.S. Powder Coated Tubular Head & Foot Bows of equal height Drop down Safety side railings on both sides Bed Mounted on PVC Stumps Pre Treated and Epoxy powder Coated Frame work of precise steel tubes Top of perforated CRCA sheet main frame of rectangular pipe with full length drop side railing Bows of equal height mounted on rubber stump Finish in epoxy powder coated approximate Size{LxWxH} {137x 76 x60} mm. | |

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003