

**Product Name :**  
Peel Off Strength Tester Machines

**Product Code :**  
NLAB-TESTINGLABM54001



**Description :**

Peel Off Strength Tester Machines

**Technical Specification :**

Peel Off Strength Tester to test the self-adhesiveness of foam sheet. The strength of adhesion of self-adhesive foam sheets to the substrate is determined by pasting a rectangular test specimen cut out of the sheet under test over a smooth glass surface, pressing it down with the help of rubber-lined metal roller, and peeling off one of its ends from the substrate at an angle of 180°; at a constant speed with the help of a motorized arrangement. & The force needed to peel off the test specimen gives a measure of the peel strength of the foam sheet under test

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003