

Product Name : Peel Strength Tester	Product Code : NLAB-TESTINGLABM54003
Description : Peel Strength Tester	
Technical Specification : Peel Strength Tester is used to determine Peel / Seal Strength of Polymer film, Plastic films, Adhesive Tapes. This machine is useful for Peel Strength, Bond Strength, Adhesion Strength, Seal Strength, Seam Strength etc. This is Table Top Model having digital Display of Load. The maximum Load (Peak value) hold on memory. This machine has self calibration option.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003