

Product Name :
Peg Board

Product Code :
NLAB-LABORATORYEM81005



Description :

Peg Board

Technical Specification :

Application: Office Building, Hospital, School, Other

Material: PP, PP

Feature: chemical resistance

Function: Air-dried test tube

Column: 27 Or 52

Size:W400*H500*D120mm W550*H700*D120mm

High-strength PP material,Corrosion resistance, acid and alkali resistance,

Drain design,Humanized design, effectively avoiding the accumulation of water pollution

27 or 52 drip column design,Detachable drip column, easier to clean

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003