

Product Name : Penta Microscope	Product Code : NLAB-EDUSCHOOLEQ1730149
Description :	
Penta Microscope	
Technical Specification :	
Condenser: Swing out condenser N.A. 09/0.25 to provide perfect illumination under all magnification to all the observers.	
Magnification: 40X-1000X Standard Optionally up to 2000X	
Eyepiece: Ultra wide fields and compensating WF 10x paired eyepiece 5 pairs.	
Observation Head: One number high Transmission and top of Trinocular head include Compensation free 360 rotatable. Four numbers compensations free Binocular Head inclined and 360 rotatable.	
Light Pointer: A Movable Built In Green Led Pointer (brightness Adjustable).	
Focusing: Low and forward position Co-axial coarse and fine focusing mechanism supported on ball bearing with tension control ring and pre-focus stopper. Coarse focus 30mm. Fine motion reading to 0.002mm.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003