

**Product Name :**  
Pentax Tribrach Adaptor

**Product Code :**  
NLABMEDICAL-M3300003

**Description :**

Pentax Tribrach Adaptor

**Technical Specification :**

Tribrach with Circular Bubble without Optical Plummet. For accepting above Adaptor.

Pentax Tribrach Adaptor with Rotating.

Optical Plummet with 90" main bubble.

60-second plate vial, adjustable.

Precision-machined axis rotates 360°.

For accepting any holder with target plate.

The over point accuracy is 0.5 mm at 1.5 m.

2.5X optical plummet.

Fits all standard three-prong tribrach.

Adjustable height:100 mm-to-120 mm.

Compatible with all prisms with a 24.9 mm bayonet base.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003