


Product Name :
Penumatic Circuit Trainer Kit (Basic)

Product Code :
NLAB-ENGINEERINGLB47004



Description :
Penumatic Circuit Trainer Kit (Basic)

Technical Specification :

SPECIFICATIONS:

COMPONENTS	
QTY.	NO
Trolley with Caster wheels	1
Quick Push-Pull connection	122

High Pressure Tubing (PU6)	5
5/2 way lever operated valve	4
5/3 way Lever operated DC valve	5
Pressure Gauge 4" Dial	6
OR Gate valve	7
Quick exhaust valve	8
One way flow control valve	9
Single acting cylinder	10
Double acting cylinder	11
FRL unit	12
Manifold	

113

SERVICE REQUIRED: - 7 Kg/cm² Air Compressor with 15amp single Phase Electric Connection

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003