

Product Name : Peptone Type III, Bacteriological	Product Code : NLAB-CHEMICALSIND20016
Description : Peptone Type III, Bacteriological	
Technical Specification : Peptone Type III, Bacteriological Peptone Type III, Bacteriological is prepared under controlled conditions by enzymatic digestion of selected fresh meat. It is recommended for general bacteriological work including mass cultivation of microorganisms. Principle And Interpretation Peptone Type III, Bacteriological is highly nutritious since, it can be used for microbiological purpose as media ingredient. It can also be used for commercial production of enzymes, vaccines, antibiotics etc. Quality Control Reaction: Reaction of 2% w/v aqueous solution at 25°C. pH: 5.0- 7.0 Solubility: Freely soluble in distilled/purified water, insoluble in alcohol. Appearance: Light yellow to brownish yellow homogenous free flowing powder characteristic odour but not putrescent.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003