

<b>Product Name :</b> Peptone Type III, Bacteriological	<b>Product Code :</b> NLAB-CHEMICALSIND20016
<b>Description :</b> Peptone Type III, Bacteriological	
<b>Technical Specification :</b> Peptone Type III, Bacteriological  Peptone Type III, Bacteriological is prepared under controlled conditions by enzymatic digestion of selected fresh meat. It is recommended for general bacteriological work including mass cultivation of microorganisms.  Principle And Interpretation  Peptone Type III, Bacteriological is highly nutritious since, it can be used for microbiological purpose as media ingredient. It can also be used for commercial production of enzymes, vaccines, antibiotics etc.  Quality Control  Reaction: Reaction of 2% w/v aqueous solution at 25°C.  pH: 5.0- 7.0  Solubility: Freely soluble in distilled/purified water, insoluble in alcohol.  Appearance: Light yellow to brownish yellow homogenous free flowing powder characteristic odour but not putrescent.	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003