

**Product Name :**  
Permeability Test System

**Product Code :**  
NLAB-VOCATIONALTVET18027



**Description :**

Permeability Test System

**Technical Specification :**

Permeability Cell:-

Cell Diameter Ø 75 mm.

The (b) type permeability cell is specially designed for constant head tests on incoherent soil samples.

It has a central body with 3 outlets for 3 manometer tubes through which water is collected after seepage through the sample.

Inlet and drainage ducts are found in the metal base.

A movable stem enables stabilization of the upper surface of the sample.

Weight 3 Kg Standpipe with three manometers tubes:-

Comprises 3 glass tubes all of the same diameter, relevant graduated rods, taps and metal pedestal for bench mounting.

Dimensions: 320 x 1450 x 500 mm.

Weight: 15 g

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003