

Product Name :
Perspirometer Testing Machines

Product Code :
NLAB-TESTINGLABM68001



Description :

Perspirometer Testing Machines

Technical Specification :

Description:

Perspirometer is used for determining the resistance of the color of textiles, fabric, yarn and garments to the action of human perspiration .

Colour Fastness against Perspiration with weight for loading 21 acrylic plates and 100 Nos. adjacent fabric.

Made of SS 136. (Fully SS)

Weight: -2nos.

Size of Specimen : 10cm x 4 cm

Force on Specimen : 5 kg/4.53 kg

Size of acrylic plates. : 12.5 cm x 6.5 cm. (As per AATCC & ISO Standards)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003