

Product Name :
Pharma Ampoule Filling Machine

Product Code :
NLAB-PHARMACYLABM12001



Description :

Pharma Ampoule Filling Machine

Technical Specification :

The machine is a new generation of the ampoule drawing ,filling and sealing combined equipment for injections .

It adopts the rectangular feeding instead of the traditional ampoule feeding methods .

In such case ,the feeding is steady with no easy broken ampoule .

A modern clutch is mouthed under the ampoule feed hopper .

Therefore ,the ampoule feeding can be stopped at any time while the machines is still running .

So the ampoule lost can be reduced .

In the part of stopping filling where changes electromagnetic switch into mechanical switch .

It is reliable and durable .

Bottle reverse and bottle arranging stricture are used in the device of bottle outlet .

It reduces the broken rate of the finished products of ampoules tremendously.

Specification: 1~2ml

Output: 4000~4800pcs/H

Power and voltage: 0.37kw , 220v/110v,

Fuel consumption: 2-2.5m³/h(gas) 0.08-0.1MPa(LPG)

Combustion-supporting Consumption: 0.7m³/h Oxygen 0.1MPa

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003