

Product Name :
Pharma Automatic Aerosol Filling Machine Production Line

Product Code :
NLAB-PHARMACYLABM13003



Description :

Pharma Automatic Aerosol Filling Machine Production Line

Technical Specification :

Application: Beverage, Chemical, Food, Medical

Automatic Grade: Automatic

Driven Type: Pneumatic

Voltage: 220V/380V (Used for air compressor)

Power: 30KW (Used for air compressor)

Dimension (L*W*H): 9000*1050*1700mm

Weight: 2500KG

Usage : Aerosol can filling

Filling capacity: 0-1000ml (can be customized)

Can diameter: 35--70mm

Applicable Aerosol Can Height: 70-330mm

Productivity: 1800-3000 cans/hr

Machine main Material: Stainless Steel 304

Applicable spray can valve: 1 inch (25.4mm)/(can be customized)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003