

**Product Name :**  
Pharma Automatic Tablet and Capsule Blister Packing Machine

**Product Code :**  
NLAB-PHARMACYLABM15002



**Description :**

Pharma Automatic Tablet and Capsule Blister Packing Machine

**Technical Specification :**

Punching frequency: 20-40times/min

Capacity: 2000 plates/ hour

Max forming area and depth: 70×100×12mm(designed according user's requirement )

Standard Plate: 80×57(designed according user's requirement)

Adjustable scope of travel: 30-80(designed according user's requirement)

Air-compressor: 0.6-0.8Mpa

Total power: 380V/220V 50Hz 1.8Kw

Main motor power: 0.75Kw

PVC: 0.15-0.5×80mm

PTP: 0.02-0.035×80mm

Dialysis paper: 50-100g×80mm

Mold cooling: Use of water cooling or chiller cooling

dimensions: 1400×460×950mm(L×W×H)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003